

Medium used for the collection, transport and preservation of clinical samples.

• CONTENTS (Liter)

Sodium Thioglycollate	1.5 g
Disodium Phosphate	1.1 g
Calcium Chloride	0.09 g
Sodium Chloride	5.0 g
Agar	5.0 g
Final pH = 8.4 \pm 0.2 at 25°C	

PROCEDURE

Suspend 12.69 G of powder in 1 L of distilled or deionized water. Heat to boiling until completely dissolved. Mix well. Stream for 15 minutes. Pour into tubes.

INTERPRETATION

Cary Blair Transport Medium is a medium used for the collection, transport and preservation of clinical samples. Sodium thioglycollate provides reducing agent. Disodium phosphate is the buffering agent. Calcium chloride supplies electrolytes for bacterial survival. Sodium chloride maintains the osmotic balance. Agar is the solidifying agent.

TECHNIC

Put the swab specimens to the medium. Transport the medium with specimen to the laboratory at room temperature as soon as possible. Inoculate the swab specimen to the proper medium. Transport and inoculation of specimens have to be performed within 48 - 72 hours. Incubate the inoculated medium at $35 \pm 2^{\circ}$ C for 18 - 48 hours under appropriate conditions. Refer appropriate references for recommended test procedure.

• QUALITY CONTROL FOR USE

 $\begin{array}{l} \underline{Dehydrated\ medium}\\ Appearance:\ free-flowing,\ homogeneous\\ Color:\ light\ beige\\ \underline{Prepared\ medium}\\ Appearance:\ clear\ to\ slightly\ opalescent\\ Color:\ very\ light\ amber\\ Incubation\ conditions:\ 35\ \pm\ 2^\circ C\ /\ 18\ -\ 48\ hours\ under\ appropriate\ conditions\ and\ medium \end{array}$

Microorganism	ATCC	Recovery
Heamophilus influenzae	19418	good
Salmonella typhimurium	14028	good
Shigella flexneri	12022	good
Streptococcus pneumoniae	6305	good
Vibrio parahaemolyticus	17802	good

STORE

The powder is very hygroscopic. Store the powder at room temperature, in a dry environment, in its original container tightly closed and use it before the expiry date on the label. Store prepared medium at 2 - 8°C.

REFERENCES

- 1. Cary, S.G., and E. B. Blair. (1964). J. Bacteriol. 88:96-98.
- 2. Cary, S.G. (1965). Am. J. Clin. Path. 43:294-296

PACKAGE

Cat. No : MB-C1403 500 G

