

# Alkaline Peptone Water (MFDS)



Medium used for the cultivation and enumeration of *Vibrio* spp. from food and water.

## • CONTENTS (Liter)

Peptone	10.0 g
Sodium Chloride	20.0 g
Final pH = 8.6 ± 0.2 at 25°C	

## • PROCEDURE

Suspend 30.0 G of powder in 1 L of distilled or deionized water. Heat to boiling until completely dissolved. Sterilize by autoclave at 121°C for 15 minutes. Cool to 45 - 50°C in water bath. Mix well. Pour into tubes.

## • INTERPRETATION

Alkaline Peptone Water is a medium used for the cultivation and enumeration of *Vibrio* spp. from food and water. The high value of pH slows down the development of common contaminants so that *Vibrio* spp. can be cultivated effectively after 6 - 7 hours of incubation. Peptone provides nitrogen and vitamin sources. Sodium chloride maintains the osmotic balance and promotes the growth of *Vibrio* spp.

## • TECHNIC

Inoculate the specimen using a sterile needle to the medium. Incubate at 36 ± 1°C for 18 - 24 hours. Refer appropriate references for recommended test procedure.

## • QUALITY CONTROL FOR USE

### Dehydrated medium

Appearance: free-flowing, homogeneous

Color: light beige

### Prepared medium

Appearance: clear to very slightly opalescent

Color: light amber

Incubation conditions: 36 ± 1°C / 18 - 24 hours

Microorganism	ATCC	Growth
<i>Vibrio parahaemolyticus</i>	17802	good

**• STORE**

The powder is very hygroscopic. Store the powder at room temperature, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration are evident. Store prepared medium at 2 - 8°C.

**• REFERENCES**

1. APHA (1985) – Standard Methods for the Examination of Water Wastewater , 16th ed.
2. Benenson, A.S., Islam M.R. & Greenough, W.B. (1964) – Rapid identification of *Vibrio cholerae* by darkfield microscopy. Bull, WHO, 30, 827.
3. Refer to the MFDS.

**• PACKAGE**

Cat. No : MB-A1098K Alkaline Peptone Water (MFDS)	500 G
--	-------